

Testing Tube Types

AST
Green
PT
Light blue
PSA
SST
CBC
Lavender
Cholesterol Select 2 answers
Red
SST
Glucose
Gray
ESR
Black
FBS
Gray
INR
Light blue
Gentamicin
Red
BC
Yellow SPS
Ddimer
Light blue
Lactic Acid
Gray
RPR
Red
BMP Please select 2 answers
Red
SST
PTT
Light blue
Estrogen
Red

Stat Lyte
Green
CMP
SST
Ca
SST
ABO/RH group
Pink
BUN
Green
Lyte
Green
Bilirubin Please select 2 answers
Red
SST
HCG Please select 2 answers
SST
Red
Tobramycin
Red
Vancomycin
Red
Alcohol
Gray
Albumin Please select 2 answers
Red
SST